



INTENDED USE :

The reagent kit is intended for the “in vitro” quantitative determination of C-Reactive Protein in serum.

CLINICAL SIGNIFICANCE:

CRP is an acute phase protein present in normal serum which increases significantly after most forms of tissue injuries, bacterial and virus infections, inflammation and malignant neoplasia. During tissue necrosis and inflammation resulting from microbial infection the CRP concentration can rise up to 300 mg/L in 12-24 hrs.

PRINCIPLE:

The CRP Turbilatex is a quantitative turbidimetric test for the measurement of C-reactive protein (CRP) in human serum or plasma. Latex particles coated with specific anti-human CRP are agglutinated when mixed with samples containing CRP. The agglutination causes an absorbance change dependent upon the CRP contents of the patient sample that can be quantified by comparison from a calibrator of known CRP concentration.

CONTENTS:

- Reagent 1: Diluent 40 ml
- Reagent 2 : Latex Antigen 10 ml
- Reagent 3 : CRP Calibrator

PREPARATION OF WORKING REAGENT

Working Reagent: Swirl the latex vial gently before use Prepare the necessary amount as follow:

- 1 ml latex Reagent + 4 ml Diluent
- CRP Calibrator : Reconstitute with 1.0 ml of Distilled water. Mix gently and keep at room temperature for 10 minutes before use.

STORAGE & STABILITY :

Working Reagent: Stable for 30 days at 2-8°C
 CRP Calibrator : Stable for 1 month at 2-8°C or 3 months at -20 C. Do not freeze; frozen latex or diluent could change the functionality of the test.

MATERIALS REQUIRED BUT NOT PROVIDED:

- Clean & Dry Glassware.
- Laboratory Glass Pipettes or Micropipettes 8rTrps.
- Bio-Chemistry Analyzer.

SAMPLES:

Fresh serum stable for 7 days at 2 -8° C or 3 months at -20° C. The samples with presence of fibrin should be centrifuged before testing. Do not use highly hemolyzed or lipemic samples.

GENERAL SYSTEM PARAMETERS

Reaction Type	Fix Time
Wave Length	546 nm (530-550 nm)
Cuvette Temperature	37°C
Delay Time	10 Sec.
Read Time	120sec
Sample Volume	5
Reagent Volume	1ml
Calibrator Concentration	AsMentioned on Calibrator vial
Zero Setting	Deionised Water
Linearity	150mg/L
Units	mg/L

PROCEDURE :

Reagent	Calibrator	Sample
Working Reagent	1ml	1ml
Calibrator	5 µl	-----
Sample	---	5 µl

CALCULATION OF RESULTS :

$$\text{CRP mg/L} = ((A2-A1)\text{Sample}/(A2-A1)\text{Calibrator}) \times \text{Calibrator Concentration}$$

NORMAL VALUE:

Up to 6 mg/L

Each laboratory should establish its own reference range.

LINEARITY:

The reaction is linear up to 150 mg/L

PERFORMANCE CHARACTERISTICS:

1. Linearity Limit: Up to 150 mg/L, under the described assay conditions.
2. Detection Limit: Values less than 2 mg/L give non-reproducible results
3. Prozone Effect: No prozone effect was detected up to 800 mg/L

QUALITY CONTROL:

It is recommended that each laboratory should prepare their own quality control scheme

LIMITATIONS AND PRECAUTIONS:

Storage conditions mentioned on the kit to be adhered.

Bilirubin (20 mg/dl) lipemia (10 g/L) and Rheumatoid factors (300 IU/ml) not interfere. Haemoglobin (>5 g/L) may interfere. Other substances may interfere

REFERENCES:

1. Lars-Olof Hanson et al. Current Opinion in Infect Diseases 1997; 10:6-201.
2. Chetana Vaishnavi. Immunology and Infectious Diseases 1996; 6: 139-144.
3. Youshitsuji Hokama et al. Journal of Clinical Lab. Status 1987; 1:15-27.
4. Kari Pulki et al. Scand J Clin Lab Invest 1986; 46: 2606-607.
5. Wemer Muller et al. Journal of Immunological Methods 1985; 80:77-90
6. Shogo Otsuji et al. Clin Chem 1982; 28/10:2121-2124.
7. Young 03. Effects of drugs on clinical laboratory test, 4th ed. AACC Press, 1995.